

Data Sheet

HUMAN NATRIURETIC PEPTIDE RECEPTOR A [NPR-A] (aa 294-308)

SYNTHETIC

Catalog No.:	AE1031AG.1 / AE1031AG.2
Name:	Human NPR-A (aa 294-308)
Synonyms:	Atrial natriuretic peptide receptor type A (ANPR-A, ANP-A), Guanylate cyclase A (GC-A)
Sequence:	LKQLKHLAYEQFNFT
Swiss-Prot No:	P16066
Gene Information:	Gene Name: NPR1, ANPRA GeneID: 4881
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	13 May 2022

For research use only

Publishing research using AE1031AG? Please let us know so that we can cite your publication as a reference.

