


## Data Sheet

# OKADAIC ACID

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	AC7380.1 / AC7380.2
<b>Immunogen:</b>	Okadaic acid coupled to a carrier peptide (FID33)
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Okadaic acid



<b>Contents:</b>	20 µl / 100 µl Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:2000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	14 April 2022

**For research use only**

**Publishing research using AC7380? Please let us know so that we can cite your publication as a reference.**

