



Data Sheet

HUMAN OSTEOCALCIN (aa 1-22)

ANTIBODY, POLYCLONAL

Catalog no.:	A 4211.1 / A 4211.2
Immunogen:	Synthetic human Osteocalcin (aa 1-22) (YLYQWLGAPVPYPDPLEPRREV)
Synonyms:	Gamma-carboxyglutamic acid-containing protein, Bone Gla protein (BGP)
Swiss-Prot No:	P02818
Gene Information:	Gene Name: BGLAP GeneID: 632
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Osteocalcin (aa 1-22) There was no cross reactivity obtained with human Osteopontin, human Osteonectin, and human Bone Sialoprotein
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:2000), ELISA (1:5000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	18 March 2022

For research use only

Publishing research using A 4211? Please let us know so that we can cite your publication as a reference.

