

Data Sheet

HUMAN OSTEOINDUCTIVE FACTOR [OIF] (aa 107-121) SYNTHETIC

Catalog no.:	AE1051AG.1 / AE1051AG.2
Name:	Human Osteoinductive Factor (aa 107-121) (RFNKIKKLTAKDFAD)
Synonyms:	Mimecan, Osteoglycin
Swiss-Prot No:	P20774
Gene Information:	Gene Name: OGN, OIF, SLRR3A GeneID: 4969
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	18 May 2022

For research use only

Publishing research using AE1051AG? Please let us know so that we can cite your publication as a reference.

