


Data Sheet

HUMAN OSTEOPONTIN

NATIVE

Catalog no.:	A 4228AG.1
Name:	Native Human Osteopontin
Synonyms:	Bone sialoprotein 1, Secreted phosphoprotein 1 (SPP-1), Urinary stone protein, Nephropontin, Uropontin
Swiss-Prot No:	P10451
Gene Information:	Gene Name: SPP1, BNSP, OPN GenelID: 6696
Source:	Human breast milk
Purification:	DEAE anion exchange and hydrophobic interaction chromatography
Storage buffer:	PBS, pH 7.3
Appearance:	Liquid
Purity:	> 90 % (SDS-PAGE)
Contents:	10 µg
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	05 May 2022



For research use only

Publishing research using A 4228AG? Please let us know so that we can cite your publication as a reference.

