

Data Sheet

HUMAN PARATHYROID HORMONE [PTH] (aa 1-34)

SYNTHETIC

Catalog No.:	A 1111AG.1 / A 1111AG.2
Name:	Human PTH (aa 1-34)
Synonyms:	Parathyrin, Parathormone
Sequence:	SVSEIQLMHNLGKHLNSMERVEWLRKKLQDVHNF
MW:	4.12 kDa
Swiss-Prot No:	P01270
Gene Information:	Gene name: PTH GeneID: 5741
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	4 May 2022

For research use only

Publishing research using A 1111AG? Please let us know so that we can cite your publication as a reference.

