

## Data Sheet

# HUMAN PARATHYROID HORMONE [PTH] (aa 1-38)

## SYNTHETIC

<b>Catalog No.:</b>	AK1105AG.1 / AK1105AG.2
<b>Name:</b>	Human PTH (aa 1-38)
<b>Synonyms:</b>	Parathyrin, Parathormone
<b>Sequence:</b>	SVSEIQLMHNLGKHLNSMERVEWLRKKLQDVHNFVALG
<b>MW:</b>	4.4 kDa
<b>Swiss-Prot No:</b>	P01270
<b>Gene Information:</b>	Gene name: PTH GenelD: 5741
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	20 May 2022

**For research use only**

**Publishing research using AK1105AG? Please let us know so that we can cite your publication as a reference.**

