

Data Sheet

RAT PARATHYROID HORMONE [PTH] (aa 53-84)

SYNTHETIC

Catalog no.:	A 1611AG.1 / A 1611AG.2
Name:	Rat PTH (aa 53-84)
Synonyms:	Parathyrin
Sequence:	KKEENVLVDGNSKSLGEGDKADV DVLVKAQSQ
MW:	3.4 kDa
Swiss-Prot No:	P04089
Gene Information:	Gene Name: Pth GenelD: 24694
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	4 May 2022

For research use only

Publishing research using A 1611AG? Please let us know so that we can cite your publication as a reference.

