

Data Sheet

HUMAN PARATHYROID HORMONE-RELATED PEPTIDE [PTHrP] (aa 53-86)

SYNTHETIC

Catalog no.:	A 1118AG.1 / A 1118AG.2
Name:	Human PTHrP (aa 53-84)
Synonyms:	Parathyroid hormone-like protein (PLP)
Sequence:	KNHPVRFSGSDDEGRYLTQETNKVETYKEQPLKTP
MW:	4.01 kDa
Swiss-Prot No:	P12272
Gene Information:	Gene Name: PTHLH, PTHRP GeneID: 5744
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	4 May 2022

For research use only

Publishing research using A 1118AG? Please let us know so that we can cite your publication as a reference.

