



Data Sheet

HUMAN PLASMINOGEN ACTIVATOR INHIBITOR 1 (PAI-1)

ANTIBODY, POLYCLONAL

Catalog no.:	A 7621.3 / A 7621.4
Immunogen:	Purified human PAI-1
Synonyms:	Endothelial plasminogen activator inhibitor, Serpin E1
Swiss-Prot No:	P05121
Gene Information:	Gene Name: SERPINE1, PAI1, PLANH1 GenelD: 5054
Host:	Goat
Matrix:	IgG, Protein G purified from serum; PBS pH 7.4
Specificity:	Human PAI-1 There was no cross-reactivity obtained with human tPA, uPA, or PAI-2.
Contents:	10 µg / 100 µg (lyophilized) Resuspend in 10 µl / 100 µl aqua bidest.
Known applications:	ELISA, RIA, immunohistochemistry (paraffin sections, cryosections) Pretreatment of paraffin sections with 0.05% Pronase E (Sigma) for 20 min. is recommended. This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	7 April 2022

For research use only

Publishing research using A 7621? Please let us know so that we can cite your publication as a reference.

