



## Data Sheet

# HUMAN PLASMINOGEN ACTIVATOR INHIBITOR 1 (PAI-1) RECOMBINANT

<b>Catalog no.:</b>	A 7621AG.1 / A 7621AG.2
<b>Name:</b>	Human PAI-1
<b>Synonyms:</b>	PAI, Endothelial plasminogen activator inhibitor, Serpin E1
<b>MW:</b>	43 kDa
<b>Swiss-Prot No:</b>	P05121
<b>Gene Information:</b>	Gene Name: SERPINE1, PAI1, PLANH1 Gene ID: 5054
<b>Source:</b>	<i>E. coli</i>
<b>Storage buffer:</b>	50 mM NaH <sub>2</sub> PO <sub>4</sub> , 100 mM NaCl, 1 mM EDTA, pH 6.6
<b>Purity:</b>	>95% by SDS-PAGE analysis
<b>Contents:</b>	10 µg / 50 µg (liquid)
<b>Store at:</b>	- 80 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	6 May 2022

**For research use only**

**Publishing research using A 7621AG? Please let us know so that we can cite your publication as a reference.**

