



## Data Sheet

# HUMAN PNEUMADIN SYNTHETIC

<b>Catalog No.:</b>	AE1134AG.1 / AE1134AG.2
<b>Name:</b>	Human Pneumadin
<b>Synonyms:</b>	PNM
<b>Sequence:</b>	AGEPKLDAGV
<b>MW:</b>	0.96 kDa
<b>Swiss-Prot No:</b>	P22103
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
<b>Store at:</b>	-20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	20 May 2022

**For research use only**

**Publishing research using AE1134AG? Please let us know so that we can cite your publication as a reference.**

