



Data Sheet

SHEEP PROLACTIN

ANTIBODY, POLYCLONAL

Catalog no.:	AS1003.1/AS1003.2
Immunogen:	Sheep Prolactin
Swiss-Prot No:	P01240
Gene Information:	Gene Name: PRL GeneID: 443317
Host:	Rabbit
Matrix:	Serum
Specificity:	Sheep, deer, and ruminants Prolactin
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	RIA (1:160 000) ¹ , immunohistochemistry This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Schams D, Rüsse I, Schallenberger E, Prokopp S, Chan JS (1984). The role of steroid hormones, prolactin and placental lactogen on mammary gland development in ewes and heifers. <i>J Endocrinol</i> 102 (1): 121-130.
Last updated on:	19 January 2023

For research use only

Publishing research using AS1003? Please let us know so that we can cite your publication as a reference.

