



Data Sheet

HUMAN RELAXIN 2 A CHAIN (aa 3-10)

ANTIBODY, POLYCLONAL

Catalog no.:	A 9220.1 / A 9220.2
Immunogen:	Synthetic human Relaxin 2 A chain (aa 3-10) (YSALANKC- β ALA-C) KLH-conjugated
Swiss-Prot No:	P04090
Gene Information:	Gene Name: RLN2 GeneID: 6019
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Relaxin 2
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	13 April 2022

For research use only

Publishing research using A 9220? Please let us know so that we can cite your publication as a reference.

