



Data Sheet

HUMAN RELAXIN 2 B CHAIN (aa 14-27)

SYNTHETIC

Catalog No.:	A 9219AG.1 / A 9219AG.2
Name:	Human Relaxin 2 B chain (aa 14-27) (ELVRAQIAICGMST- β Ala-C), sequence equivalent to human Relaxin 1 B chain (aa 16-29)
Swiss-Prot No:	P04090
Gene Information:	Gene Name: RLN2 GeneID: 6019
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	-20°C Repeated thawing and freezing must be avoided
Last updated on:	06 May 2022

For research use only

Publishing research using A 9219AG? Please let us know so that we can cite your publication as a reference.

