



## Data Sheet

# HUMAN RESISTIN (aa 33-90)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 8029.1 / A 8029.2
<b>Immunogen:</b>	Synthetic human Resistin (aa 33-90) bTG- conjugated (CQSVTSRGDLATCPRGFAVTGCTCGSACGSWDVRAETTCHCQCAGMDWTGAR CCRVQP)
<b>Synonyms:</b>	Adipose tissue-specific secretory factor, ADSF, C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2, Cysteine-rich secreted protein FIZZ3
<b>Swiss-Prot No:</b>	Q9HD89
<b>Gene Information:</b>	Gene Name: RETN, FIZZ3, HXCP1, RSTN GenelD: 56729
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human Resistin
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	13 April 2022

### For research use only

**Publishing research using A 8029? Please let us know so that we can cite your publication as a reference.**

