

Data Sheet

HUMAN RESISTIN (aa 33-90)

SYNTHETIC

Catalog no.:	A 8029AG.1 / A 8029AG.2
Name:	Human Resistin (aa 33-90)
Synonyms:	Adipose tissue-specific secretory factor, ADSF, C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein, Cysteine-rich secreted protein A12-alpha-like 2, Cysteine-rich secreted protein FIZZ3
Sequence:	<p>CQSVTSRGLATCPRGFAVTGCTCGSACGSWDVRAETTCHCQCAGMDWTGAR</p> <p>CCRVQP</p>
Swiss-Prot No:	Q9HD89
Gene Information:	<p>Gene Name: RETN, FIZZ3, HXCP1, RSTN</p> <p>Gene ID: 56729</p>
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	<p>10 µg / 100 µg (lyophilized)</p> <p>Resuspend in aqua bidest for required concentration.</p>
Store at:	<p>- 20 °C</p> <p>Repeated thawing and freezing must be avoided</p>
Last updated on:	6 May 2022

For research use only

Publishing research using A 8029AG? Please let us know so that we can cite your publication as a reference.

