



## Data Sheet

# RETINOL-BINDING PROTEIN 4 (RBP4) (aa 1-183)

## ANTIBODY, POLYCLONAL

<b>Catalog no.:</b>	A 1059.1 / A 1059.2
<b>Immunogen:</b>	Native human RBP4 (aa 1-183) full length (signal peptid cleaved)
<b>Synonyms:</b>	Plasma retinol-binding protein (PRBP)
<b>Swiss-Prot No:</b>	P02753
<b>Gene Information:</b>	Gene Name: RBP4 Gene ID: 5950
<b>Host:</b>	Rabbit
<b>Matrix:</b>	Serum
<b>Specificity:</b>	Human RBP4 (aa 1-183)
<b>Contents:</b>	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
<b>Known applications:</b>	ELISA (1:2000)  This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
<b>Store at:</b>	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

**Last updated on:** 16 March 2022

### For research use only

**Publishing research using A 1059? Please let us know so that we can cite your publication as a reference.**

