



Data Sheet

RETINOL-BINDING PROTEIN 4 (RBP4) (aa 169-176)

ANTIBODY, POLYCLONAL

Catalog no.:	A 1058.1 / A 1058.2
Immunogen:	Synthetic human RBP4 (aa 169-176) KLH-conjugated (VHNGYCDG)
Synonyms:	Plasma retinol-binding protein (PRBP)
Swiss-Prot No:	P02753
Gene Information:	Gene Name: RBP4 Gene ID: 5950
Host:	Rabbit
Matrix:	Serum
Specificity:	Human RBP4 (aa 169-176)
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA (1:2000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided

Last updated on: 16 March 2022

For research use only

Publishing research using A 1058? Please let us know so that we can cite your publication as a reference.

