

## Data Sheet

# HUMAN PROTHYMOSIN $\alpha$ (aa 75-109)

## SYNTHETIC

<b>Catalog no.:</b>	A 9540AG.1 / A 9540AG.2
<b>Name:</b>	Human ProThymosin $\alpha$ (aa 75-109)
<b>Sequence:</b>	GDEDEEAESATGKRAAEDDEDDVDTKKQKTDEDD
<b>MW:</b>	12.05 kDa
<b>Swiss-Prot No:</b>	P06454
<b>Gene Information:</b>	Gene Name: PTMA, TMSA GeneID: 5757
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	10 $\mu$ g/100 $\mu$ g (lyophilized) Resuspend in aqua bidest for required concentration
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	11 May 2022

**For research use only**

**Publishing research using A 9540AG? Please let us know so that we can cite your publication as a reference.**

