



Data Sheet

HUMAN THYMOSIN β 10 (aa 1-12)

ANTIBODY, POLYCLONAL

Catalog no.:	A 9562.1 / A 9562.2
Immunogen:	Synthetic human Thymosin β 10 (aa 1-12) BSA conjugated (ACADKPDMGEIASF)
Swiss-Prot No:	P63313
Gene Information:	Gene Name: TMSB10, PTMB10, THYB10 GenelD: 9168
Host:	Rabbit
Matrix:	Serum
Specificity:	Synthetic human Thymosin β 10 (aa 1-12), human Thymosin β 10 There was no cross reactivity obtained with human Thymosin β 4, Thymosin β 9, Thymosin β 15, Thymosin β 4 peptide (aa 1-4), and Thymosin β 4 (aa 1-14).
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA (1:1000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Hörger S, Gallert B, Echner H, Voelter W (1992). Synthese eines Thymosin beta-10-Fragments zur Entwicklung spezifischer Antikörper. <i>Z Naturforsch</i> 47b :1170-4.
Last updated on:	14 April 2022

For research use only

Publishing research using A 9562? Please let us know so that we can cite your publication as a reference.

