



Data Sheet

HUMAN THYMOSIN β 15 (aa 1-45)

ANTIBODY, POLYCLONAL

| | |
|----------------------------|--|
| Catalog no.: | A 9575.1 / A 9575.2 |
| Immunogen: | Synthetic human Thymosin β 15 (aa 1-45) KLH conjugated (MGDRPDLSEVERFDKSKLKKTTITEVKNTLPSKETIEQEKEFKRS) |
| Synonyms: | NB thymosin beta, Thymosin-like protein 8 |
| Swiss-Prot No: | P0CG34 |
| Gene Information: | Gene Name: TMSB15A, TMSL8, TMSNB GeneID: 11013, 286527 |
| Host: | Rabbit |
| Matrix: | Serum |
| Specificity: | Human Thymosin β 15 (aa 1-45) |
| Contents: | 20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest. |
| Known applications: | ELISA (1:80 000) This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls. |
| Store at: | 2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided |
| Last updated on: | 14 April 2022 |

For research use only

Publishing research using A 9575? Please let us know so that we can cite your publication as a reference.

