



Data Sheet

HUMAN THYMOSIN β 4 (aa 1-4)

ANTIBODY, POLYCLONAL

Catalog no.:	A 9521.1 / A 9521.2
Immunogen:	Synthetic human Thymosin β 4 (aa 1-4) KLH conjugated (AcSDKP)
Synonyms:	T beta-4, Fx
Swiss-Prot No:	P62328
Gene Information:	Gene Name: TMSB4X, TB4X, THYB4, TMSB4 GeneID: 7114
Host:	Rabbit
Matrix:	Serum
Specificity:	Human Thymosin β 4 (aa 1-4) There was no cross reactivity obtained with human Thymosin Tb9, Tb9 (aa 1-14), ADKP, AcGDKP, AcTDKP, AcSPKD, AcVDKP, BSA
Contents:	20 μ l / 100 μ l (lyophilized) Resuspend in 20 μ l / 100 μ l aqua bidest.
Known applications:	ELISA (1:4000) ¹ This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); -20 °C (dissolved) Repeated thawing and freezing must be avoided
References:	1. Becker S, Armbruster FP, Muller B, Echner H, Kapurnotu A, Livaniou E, Mihelic M, Stoeva S, Voelter W (1994). Theoretical and experimental epitope mapping of thymosin beta 4. <i>J Immunol Methods</i> 177(1-2): 131-137.
Last updated on:	13 April 2022

For research use only

Publishing research using A 9521? Please let us know so that we can cite your publication as a reference.

