

## Data Sheet

# HUMAN THYMOSIN $\beta_4$ SULFOXIDE (C-TERMINAL BIOTINYLATED) SYNTHETIC

<b>Catalog No.:</b>	A 9853AG.1 / A 9853AG.2
<b>Name:</b>	Human Thymosin $\beta_4$ sulfoxide (C-terminal biotinylated)
<b>Synonyms:</b>	T beta-4, Fx
<b>Sequence:</b>	SDKPD-M(O)-AEIEKFDKSKLKKTTETQEKNPLPSKETIEQEKGAGES- <u>L</u> -Biotin Underlined amino acids were added or modified for technical reasons.
<b>MW:</b>	5.2 kDa
<b>Swiss-Prot No:</b>	P62328
<b>Gene Information:</b>	Gene Name: TMSB4X, TB4X, THYB4, TMSB4 GenelD: 7114
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 $\mu$ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
<b>Store at:</b>	- 20°C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	11 May 2022

**For research use only**

**Publishing research using A 9853AG? Please let us know so that we can cite your publication as a reference.**

