

Data Sheet

HUMAN THYMOSIN β_4 SULFOXIDE (N-TERMINAL BIOTINYLATED) SYNTHETIC

Catalog No.:	A 9852AG.1 / A 9852AG.2
Name:	Human Thymosin β_4 sulfoxide (N-terminal biotinylated)
Synonyms:	T beta-4, Fx
Sequence:	Biotin-SDKPD-M(O)-AEIEKFDKSKLKKTTETQEKNPLPSKETIEQEKQAGES
MW:	5.1 kDa
Swiss-Prot No:	P62328
Gene Information:	Gene Name: TMSB4X, TB4X, THYB4, TMSB4 GeneID: 7114
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 μ g / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20°C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9852AG? Please let us know so that we can cite your publication as a reference.

