



Data Sheet

HUMAN THYMULIN

ANTIBODY, POLYCLONAL

Catalog no.: A 9530.1 / A 9530.2

Immunogen: Synthetic human Thymulin (EAKSQGGSN)

Synonyms: Thymic Factor

Host: Rabbit

Matrix: Serum

Specificity: Human Thymulin

There was no cross reactivity obtained with human Thymosin α_1 and human Thymosin β_4

Contents: 20 μ l / 100 μ l (lyophilized)

Resuspend in 20 μ l / 100 μ l aqua bidest.

Known applications: RIA (1:800), ELISA (1:15 000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)

Repeated thawing and freezing must be avoided

Last updated on: 26 January 2024

For research use only

Publishing research using A 9530? Please let us know so that we can cite your publication as a reference.

