



Data Sheet

HUMAN THYMULIN SYNTHETIC

Catalog no.:	A 9530AG.1 / A 9530AG.2
Name:	Human Thymulin
Synonyms:	Thymic Factor
Sequence:	EAKSQGGSN
MW:	0.88 kDa
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	11 May 2022

For research use only

Publishing research using A 9530AG? Please let us know so that we can cite your publication as a reference.

