

## Data Sheet

# HUMAN TRANSCRIPTION FACTOR SP7 (aa 220-234) SYNTHETIC

<b>Catalog no.:</b>	A 4213AG.1 / A 4213AG.2
<b>Name:</b>	Synthetic human Transcription Factor SP7 (aa 220-234)
<b>Synonyms:</b>	Zinc finger protein osterix
<b>Sequence:</b>	PQHVLPQDVYKPKAV
<b>MW:</b>	1.72 kDa
<b>Swiss-Prot No:</b>	Q8TDD2
<b>Gene Information:</b>	Gene Name: SP7 Gene ID: 121340
<b>Source:</b>	Synthetic
<b>Purity:</b>	> 95% (HPLC)
<b>Contents:</b>	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
<b>Store at:</b>	- 20 °C Repeated thawing and freezing must be avoided
<b>Last updated on:</b>	4 May 2022

**For research use only**

**Publishing research using A 4213AG? Please let us know so that we can cite your publication as a reference.**

