

Data Sheet

HUMAN TRANSFORMING GROWTH FACTOR β 1 (TGF β 1)

RECOMBINANT

| | |
|--------------------------|---|
| Catalog no.: | A 9600AG.1 / A 9600AG.2 |
| Name: | Human TGF β 1 (aa 250-361) (ALDTNYCFSSTEKNCCVRQLYIDFRKDLGWKWIHEPKGYHANFCLGPCPYIWSLDTQYSKVLALYNQHNP GASAAPCCVPQALEPLPIVYVYVGRKPKVEQLSNMIVRSC KCS) |
| Swiss-Prot No: | P01137 |
| Gene Information: | Gene Name: TGF β 1, TGF β GenelD: 7040 |
| Source: | HEK293 cells |
| Storage buffer: | 0.1% trifluoroacetic acid |
| Purity: | > 98% (SDS-PAGE and HPLC) |
| Contents: | 2 μ g / 10 μ g (lyophilized) Resuspend in 10 mM citric acid pH 3.0 to a concentration of 0.1-1.0 mg/ml. Do not vortex. For extended storage dilute further in a buffer containing a carrier protein (for example 0.1% BSA). |
| Store at: | - 20 °C Repeated thawing and freezing must be avoided |
| Last updated on: | 11 May 2022 |

For research use only

Publishing research using A 9600AG? Please let us know so that we can cite your publication as a reference.

