



Data Sheet

HUMAN TRIIODOTHYRONINE (T3)

ANTIBODY, POLYCLONAL

Catalog no.: A1036.1 / A1036.2

Immunogen: Synthetic T3 BSA-conjugated

Host: Rabbit

Matrix: Serum

Specificity: T3

There was no cross reactivity detected with T4.

Contents: 20 µl / 100 µl (lyophilized)

Resuspend in 20 µl / 100 µl aqua bidest.

Known applications: RIA (1:1000)

This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.

Store at: 2-8 °C (lyophilized); - 20 °C (dissolved)

Repeated thawing and freezing must be avoided

Last updated on: 16 March 2022

For research use only

Publishing research using A 1036? Please let us know so that we can cite your publication as a reference.

