

Data Sheet

HUMAN URODILATIN / CDD-ANP (aa 95-104)

SYNTHETIC

Catalog no.:	A 8021AG.1 / A 8021AG.2
Name:	Human Urodilatin (aa 95-104) (TAPRSLRRSS)
Synonyms:	CDD-ANP, PREPRONATRIODILATIN
Swiss-Prot No:	P01160
Gene Information:	Gene Name: NPPA, ANP, PND GeneID: 4878
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (10 µg/µl) liquid Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	06 May 2022

For research use only

Publishing research using A 8021AG? Please let us know so that we can cite your publication as a reference.