

Data Sheet

RAT URODILATIN / CDD-ANP (aa 95-126)

SYNTHETIC

Catalog no.:	A 8023AG.1 / A 8023AG.2
Name:	Rat Urodilatin / CDD-ANP (aa 95-126)
Synonyms:	CDD-ANP, Prepronatriodilatin
Sequence:	AGPRSLRALGRGPWDPSDRSALLKSKLRALLA
MW:	3.43 kDa
Swiss-Prot No:	P01161
Gene Information:	Gene Name: Nppa GeneID: 24602
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration
Store at:	- 20 °C Repeated thawing and freezing must be avoided
References:	1. Kirchoff K, Bub A, Marxen P, Forssmann WG (1994). Urodilatin inhibits sodium reabsorption in the isolated perfused rat kidney. <i>Exp Nephrol</i> 2(6): 351-357. 2. Kirchoff K, Forssmann WG (1995). Release of urodilatin from perfused rat kidney and from cultured neonatal rat kidney cells. <i>Pflugers Arch</i> 430(5): 739-744.
Last updated on:	6 May 2022

For research use only

Publishing research using A 8023AG? Please let us know so that we can cite your publication as a reference.

