



Data Sheet

HUMAN UROKINASE PLASMINOGEN ACTIVATOR [uPA]

ANTIBODY, POLYCLONAL

Catalog no.:	A 7611.1 / A 7611.2
Immunogen:	Purified human urokinase plasminogen activator [uPA]
Synonyms:	Urokinase-type plasminogen activator, EC=3.4.21.73
Swiss-Prot No:	P00749
Gene Information:	Gene Name: PLAU GeneID: 5328
Host:	Goat
Matrix:	Serum
Specificity:	Human Urokinase plasminogen activator There was no cross-reactivity obtained with human tPA or human PAI-2.
Contents:	20 µl / 100 µl (lyophilized) Resuspend in 20 µl / 100 µl aqua bidest.
Known applications:	ELISA, immunohistochemistry (paraffin sections, 1:50; cryosections, 1:50-1:100) Recommended procedure for paraffin sections: Pretreatment with 0.05% Pronase E (Sigma) for 20 minutes. This antibody has not been tested for use in all applications. This does not necessarily exclude its use for non-tested procedures. The stated dilutions are recommendations only. We suggest that the applicant titrates the antibody in his/her system using appropriate negative/positive controls.
Store at:	2-8 °C (lyophilized); - 20 °C (dissolved) Repeated thawing and freezing must be avoided
Last updated on:	7 April 2022

For research use only

Publishing research using A 7611? Please let us know so that we can cite your publication as a reference.

