

Data Sheet

HUMAN VASOACTIVE INTESTINAL PEPTIDE-RECEPTOR 2 [VIP-R2]

SYNTHETIC

Catalog no.:	AE6302AG.1 / AE6302AG.2
Name:	Human VIP-R2 (aa 3-14)
Synonyms:	Helodermin-preferring VIP receptor, Pituitary adenylate cyclase-activating polypeptide type III receptor (PACAP-R3), VPAC2
Sequence:	RFHLEIQEEETK
MW:	1.56 kDa
Swiss-Prot No:	P41587
Gene Information:	Gene Name: VIPR2, VIP2R Gene ID: 7434
Source:	Synthetic
Purity:	> 95% (HPLC)
Contents:	100 µg / 1 mg (lyophilized) Resuspend in 0.1% acetic acid for required concentration.
Store at:	- 20 °C Repeated thawing and freezing must be avoided
Last updated on:	19 May 2022

For research use only

Publishing research using AE6203AG? Please let us know so that we can cite your publication as a reference.

