

Data Sheet

HUMAN VITAMIN D BINDING PROTEIN

ANTIBODY, POLYCLONAL

Catalog no.: A 1020.3 / A 1020.4

Immunogen: Human GC Globulin

Synonyms: Vitamin D-binding protein (VDBP, DBP), Gc protein-derived macrophage activating factor1 (Gc-MAF, GcMAF1); Gc-globulin, Group-specific component (Gc1), Vitamin D-binding protein-macrophage activating factor (DBP-maf)

Swiss-Prot No: P02774

Gene Information: Gene Name: GC
GeneID: 2638

Host: Rabbit

Matrix: IgG; purified from serum, 100mM NaCl, 15 mM NaN₃

Specificity: GC Globulin

Contents: 10 µg / 100 µg (liquid)

Known applications: ELISA

This antibody has not been tested for use in all applications. This does not necessarily exclude its use in non-tested procedures. The stated dilutions are recommendations only. End users should determine optimal dilutions in their system using appropriate negative/positive controls.

Store at: - 20 °C

Last updated on: 04 April 2025

For research use only

Publishing research using A 1020? Please let us know so that we can cite your publication as a reference.

